

Tech Tips that you can use today!

To ensure the best case possible for you and your patient we wanted to let you know the following..

Impression working time and total set time

1. Always check the working and set time of your light body and the heavy body impression materials.
2. If the light body material is one, (1) minute or less it is not recommended to use for more than one unit and your timing is very important between light body and the heavy body. If the timing is off it can cause a remake.
3. Working time of the light body of two, (2) minutes or more is recommended for two to four units. Again the timing between the light body and the heavy body is important to prevent remakes.
4. The total set time must be followed. If it is (4) minutes leave the impression in for (5) minutes. This is insurance.
5. Never leave the patient alone with the impression in the mouth. Always be there to lightly hold it in place.
6. Always use impression tray adhesive.